slit_conrner_detection

June 19, 2019

1 Overview


This notebook is used for the development of the auto slit conrner detection algorithm, which is critical for the automated detector drift correction.

2 Initial test with the first image

let's import one images and start with the code we submitted to tomopy a long time ago.


```
[3]: rawimg = plt.imread('data/test_midregion_1.tif')
fig, axes = plt.subplots(1,2,figsize=(10, 5))
img = exposure.equalize_hist(rawimg)
axes[0].imshow(rawimg, 'jet')
axes[1].imshow(img, 'jet')
```

[3]: <matplotlib.image.AxesImage at 0x1c19624748>



```
[4]: tmp = np.log(medfilt2d(img.astype(float))+1)
    col_prof = gaussian_filter1d(np.average(tmp, axis=0), sigma=11)
    dot_col_prof = np.gradient(col_prof)
    left = np.argmax(dot_col_prof)
    right = np.argmin(dot_col_prof)
    row_prof = gaussian_filter1d(np.average(tmp, axis=1), sigma=11)
    dot_row_prof = np.gradient(row_prof)
    top = np.argmax(dot_row_prof)
    bot = np.argmin(dot_row_prof)
    fig, axes = plt.subplots(1, 2, figsize=(15, 10))
    axes[0].plot(dot_col_prof)
    axes[0].plot(left, dot_col_prof[left], 'ro')
    axes[0].plot(right, dot_col_prof[right], 'ro')
    axes[1].plot(dot_row_prof)
    axes[1].plot(top, dot_row_prof[top], 'ro')
    axes[1].plot(bot, dot_row_prof[bot], 'ro')
```


[4]: [<matplotlib.lines.Line2D at 0x1c1a03af28>]



```
[5]: fig, ax = plt.subplots(1,1,figsize=(18, 10))
    ax.imshow(rawimg, 'jet')
    ax.axvline(left, color='m')
    ax.axvline(right, color='r')

ax.axvline(top, color='g')
    ax.axhline(bot, color='w')
```

[5]: <matplotlib.lines.Line2D at 0x1c1a68af60>



Making it into a function

```
[7]: def guess_slit_box(img):
        # Contrast stretching
       pl, ph = np.percentile(img, (2, 98))
       img = exposure.rescale_intensity(img, in_range=(pl, ph))
        # equilize hist
        img = exposure.equalize_adapthist(img)
        # map to log to reveal transition box
       img = np.log(medfilt2d(img.astype(float))+1)
        # get row and col profile gradient
       pdot_col = np.gradient(gaussian_filter1d(np.average(img, axis=0), sigma=11))
       pdot_row = np.gradient(gaussian_filter1d(np.average(img, axis=1), sigma=11))
       return {
            'left': np.argmax(pdot_col),
            'right': np.argmin(pdot_col),
            'top': np.argmax(pdot_row),
            'bot': np.argmin(pdot_row),
       }
```

```
# ------
# testing
rawing = plt.imread('data/test_midregion_1.tif')

edges = guess_slit_box(rawing)


fig, ax = plt.subplots(1,1,figsize=(18, 10))

# ax.imshow(rawing, 'gray')
# ax.imshow(exposure.equalize_hist(rawing), 'gray')
ax.imshow(exposure.equalize_adapthist(rawing), 'gray')

ax.axvline(edges['left'], color='r')
ax.axvline(edges['right'], color='m')
ax.axvline(edges['top'], color='g')
ax.axhline(edges['bot'], color='b')

print(edges)
```

{'left': 461, 'right': 1335, 'top': 350, 'bot': 1069}



3 Systematic test under different conditions


3.1 Standard case

In most cases, the majority of detector should be used for imaging.

```
fig, axs = plt.subplots(4, 1, figsize=(15, 30))

for i in range(4):
    img = plt.imread(f'data/test_bigregion_{i}.tif')
    edges = guess_slit_box(img)
    axs[i].imshow(exposure.equalize_adapthist(img), 'gray')
    axs[i].axvline(edges['left'], color='r')
    axs[i].axvline(edges['right'], color='m')
    axs[i].axhline(edges['top'], color='g')
    axs[i].axhline(edges['bot'], color='b')

plt.subplots_adjust(wspace=0, hspace=0.05)
```




3.2 Mid region case

In some cases, a good region of the detector is purposely selected for imaging.

```
fig, axs = plt.subplots(4, 1, figsize=(15, 30))

for i in range(4):
    img = plt.imread(f'data/test_midregion_{i}.tif')
    edges = guess_slit_box(img)
    axs[i].imshow(exposure.equalize_adapthist(img), 'gray')
    axs[i].axvline(edges['left'], color='r')
    axs[i].axvline(edges['right'], color='m')
    axs[i].axhline(edges['top'], color='g')
    axs[i].axhline(edges['bot'], color='b')

plt.subplots_adjust(wspace=0, hspace=0.05)
```




3.3 Small region cases

In very rare situation, the imaging region is shrinked to a really small region on the detector.

```
fig, axs = plt.subplots(4, 1, figsize=(15, 30))

for i in range(4):
    img = plt.imread(f'data/test_smallregion_{i}.tif')
    edges = guess_slit_box(img)
    axs[i].imshow(exposure.equalize_adapthist(img), 'gray')
    axs[i].axvline(edges['left'], color='r')
    axs[i].axvline(edges['right'], color='m')
    axs[i].axhline(edges['top'], color='g')
    axs[i].axhline(edges['bot'], color='b')

plt.subplots_adjust(wspace=0, hspace=0.05)
```




Basically, as long as all four blade slits are used, the function should be able to return the approximated location of the four blades without any trouble. However, if there is no slits at all, random results will be returned.

```
fig, axs = plt.subplots(4, 1, figsize=(15, 30))

for i in range(4):
    img = plt.imread(f'data/test_noslit_{i}.tif')
    edges = guess_slit_box(img)
    axs[i].imshow(exposure.equalize_adapthist(img), 'gray')
    axs[i].axvline(edges['left'], color='r')
    axs[i].axvline(edges['right'], color='m')
    axs[i].axhline(edges['top'], color='g')
    axs[i].axhline(edges['bot'], color='b')

plt.subplots_adjust(wspace=0, hspace=0.05)
```



The key point here is that we should only use this slit box finder if we know that slits are used
in the experiment, not the other way around.